



ADJUSTABLE FALL MANAGEMENT BED RAIL

Model #: EZBR-2 • Made in China



ASSEMBLY INSTRUCTIONS

CAUTION: Please inspect all parts carefully for cracks/breaks prior to use.

CAUTION: Please inspect all parts carefully for cracks/breaks prior to use.

1. Insert handle into cross bars to desired height; be sure adjustment buttons pop out. Be sure buttons are set evenly on both sides.
2. Insert bottom stability bar into cross bars simultaneously until adjustment button pops out.
3. Maximum weight limit: 300 lbs.

CAUTION: Be sure push buttons are aligned equally. Push buttons must fully protrude through the adjustment hole to ensure a positive lock.

CAUTION: POSSIBLE INJURY. Make sure side rails are secured properly before using bed.

CAUTION: RISK OF SERIOUS INJURY. Individuals with physical limitations who cannot prevent themselves from rolling or climbing out of the bed may require other means of safe positioning.

CAUTION: Side rails do not fall within any weight limitations and may be damaged if excessive pressure is placed on them. Side rails are for the purpose of aiding in turning and repositioning within the bed and are not intended to be used for restraint purposes.

STRAP ASSEMBLY

Safety strap provides extra stability, especially for hospital and twin beds where the mattresses are lighter and smaller. 84" safety strap will fit up to King size beds. 1. Insert bed rail under mattress.

2. Go to the opposite side of the bed and lift mattress until bottom stability bar of bed rail is visible.
3. Place strap around bed rail bottom stability bar and around bed frame bar directly below. Insert buckle ends of strap and adjust until strap is secured with a snug fit.

USER INSTRUCTIONS

Side rails are for the purpose of aiding in turning and repositioning within the bed and are not intended to be used for restraint purposes.

This product can be used as an aid in balance and support when getting in and out of the bed. This product is NOT designed to support the total weight of the user.

CARE INSTRUCTIONS

Periodically check that all push buttons are securely engaged and strap is secure.

One Year Limited Warranty

Personal Safety Corporation warrants this product to be free from factory defects and workmanship for a period of one year from date of purchase. Void if this product is misused, abused altered, modified or used for purposes other than as an adult bed rail. No warranty of use expressed or implied.

CAUTION! This product is not a restraint. It will not prevent a resident from getting out of their bed. This product is to be used in conjunction with a total fall prevention program.



Please see our complete line of Secure[®] fall/wandering prevention & patient safety products at www.SecureSafetySolutions.com

Highest Quality...Lowest Prices...Guaranteed!

Personal Safety Corporation • Hiawatha, Iowa 52233 • 800.373.2873