



# ELESS OTION SENSOR

Model #: CAS-PIR • Made in PRC

# PARTS AND FUNCTION





#### BATTERY INSTALLATION






**PAIRING** 

The CAS-PIR motion sensor must be paired with the CAS-PGR Wireless Caregiver Pager and/or CAS-WCR Wireless Caregiver Receiver before using. Always test the system before each use.

- Step 1: Insert battery and switch on the PIR sensor (Figure 1).
- Step 2: Insert batteries into the CAS-PGR wireless pager (Figure 2) and/or CAS-WCR wireless receiver (Figure 3).



Zone selecto


Melody

Figure 3

button

- Step 3: Enter pager and/or receiver pairing mode by pressing and holding the zone select button • for a few seconds until zone light "I" begins to flash.
- Step 4: Press the zone select button again to switch to the zone you wish to pair with. Press melody selector 🞵 to scroll through the melodies. Once you hear your preferred

melody, do not press again as it is now selected.




### TEST THE SENSOR

Insert the battery and switch on the PIR sensor, it will start to detect motion after 30 seconds. Walk to the front of the sensor within 26' 2-61/64 (8m).

Once motion has been detected.the sensor will send a signal to the CAS-PGR pager and/ or CAS-WCR receiver, then go into sleep mode for 5 seconds or 5 minutes, depending on the timer settina.



# PIR SENSOR EFFECTIVE AREA



# ONE YEAR LIMITED WARRANTY

Personal Safety Corporation warrants this product to be free from factory defects and workmanship for a period of one year from date of purchase. Warranty is void if this product is misused, abused or damaged. No warranty or responsibility of use expressed or implied.

FCC STATEMENT: THIS EQUIPMENT HAS BEEN TESTED AND FOUND TO COMPLY WHITTHE LIMITS FOR A CLASS B DIGITAL DEVICE, PURSUANT TO PART 15 OF THE FCC RULES. THESE LIMITS ARE DESIGNED TO PROVIDE REASONABLE PROTECTION AGAINST HARMFUL INTERFERENCE IN A RESIDENTIAL INSTALLATION. THIS EQUIPMENT GENERATES, USES AND CAN RADIATE RADIO FREQUENCY ENERGY AND, IF NOT INSTALLED AND USED IN ACCORDANCE WITH THE INSTRUCTIONS, MAY CAUSE HARMFUL INTERFERENCE TO RADIO COMMUNICATIONS.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- · Reorient or relocate the receiving antenna.
- · Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation

Note: Modifications to this product will void the user's authority to operate this equipment.



Please see our complete line of Secure® fall/wandering prevention & patient safety products at www.SecureSafetySolutions.com

Highest Quality...Lowest Prices...Guaranteed!

Personal Safety Corporation • Hiawatha, Iowa 52233 • 800.373.2873