

# WHEELCHAIR OR WALKER POUCH

Model # WWP-1B • Made in PRC • One Year Warranty

The Secure<sup>®</sup> Wheelchair or Walker Pouch is designed to hang on the inside of a wheelchair armrest, beneath the wheelchair seat, or on a standard walker to keep essential items close at hand. The pouch may also be hung on the outside of most power chair armrests. Easily and securely attaches using Velcro fasteners.

## Product Features:

- Main zippered pocket holds essential personal items such as a wallet, brush/comb, and other small items
- Fleece lined pocket is ideal for cell phones, eyeglasses, sunglasses or other small items
- Expanding front pocket keeps other necessities close at hand
- Easily attaches to any wheelchair or walker
- One year warranty

## Installation:

To install the Secure<sup>®</sup> Wheelchair or Walker Pouch, simply unhook both Velcro fastening straps and wrap each fastening strap around a wheelchair armrest bar or front/side of walker frame (see fig. 1). Make sure the Velcro fasteners do not interfere with the armrest surface. The pouch may also be attached by wrapping the fasteners to the wheelchair frame underneath the seat.




Figure 1

**Personal Safety Corporation warrants this product to be free from factory defects and workmanship for a period of one year from date of purchase. Void if product's outer covering is misused, abused, cut or punctured. No warranty or responsibility to use expressed or implied.**

Visit [www.SecureSafetySolutions.com](http://www.SecureSafetySolutions.com) for more  
Wheelchair Cushions, Belts & Accessories



Secure<sup>®</sup> Convex Pommel Cushion



Secure<sup>®</sup> Gel Cushion



Secure<sup>®</sup> Deluxe Arm Support



Please see our complete line of Secure<sup>®</sup> fall/wandering prevention & patient safety products  
at [www.SecureSafetySolutions.com](http://www.SecureSafetySolutions.com)

**Highest Quality...Lowest Prices...Guaranteed!**

**Personal Safety Corporation • Hiawatha, Iowa 52233 • 800.373.2873**